

# THE CLINISPIN™ CT20

## A COMPACT, VARIABLE SPEED CENTRIFUGE



- High speed and g-force
- Complete with rotor for 6 x 10/15ml tubes
- Adapters for common blood draw tubes available
- Easy to use
- Variable speed in 100 RPM increments
- Compact design
- Lightweight

### Suggested Applications

- Doctors offices
- Small laboratories
- University departments
- Coagulation surgeries
- Contract research / pharmaceutical organisations
- Veterinary surgeries

### Specifications

Speed Control	Digital - variable (100 RPM increments)
<b>Max. Speed</b>	<b>6,500 RPM</b>
Max RCF	4,000 x g
<b>Max. Capacity</b>	<b>6 x 15ml</b>
Timer	30 sec to 30 min with continuous hold features
<b>Acceleration Time</b>	<b>14 seconds</b>
Deceleration Time	27 seconds
<b>Rotor Angle</b>	<b>30°</b>
Voltage	120V or 230V
<b>Frequency</b>	<b>50 - 60 Hz</b>
Weight	4 kg
<b>Dimensions</b>	<b>18h x 21w x 24d cm</b>

Visit [www.woodleylabdiagnostics.com](http://www.woodleylabdiagnostics.com) for full product range