

Clinispin horizon 6 Flex FA & 6 FA

Fixed Angle Routine Centrifuges

Spin 125mm tubes and meet fixed angle processing needs for PPP, PRP, urine and more with the Clinispin horizon Fixed Angle Routine Centrifuges



- Fixed Angle Rotor
- Clear Shatterproof Lid for Easy Monitoring
- Lid Safety Interlock System
- Brushless DC Motor
- Cool Operation

Clinispin horizon 6 Flex FA

Digital display and 10 programmable cycles

Clinispin horizon 6 FA

Simple two-button operation to select a preset cycle or lock to a single cycle

Suggested Applications

- Doctor's Offices
- Small Laboratories
- Veterinary Surgeries
- Emergency Rooms
- Contract Research/
Pharmaceutical Organisations

Visit www.woodleylabdiagnostics.com for full product range

Clinispin horizon 6 Flex FA & 6 FA

Fixed Angle Routine Centrifuges

Spin 125mm tubes and meet fixed angle processing needs for PPP, PRP, urine and more with the Clinispin horizon Fixed Angle Routine Centrifuges



Lid lights indicate cycle status:
ready, running, complete

LED light indicates current cycle

View cycle and remaining spin
time on screen

Program your validated cycle and
give your cycle a custom name

Store up to 10 programs or use
default presets for chemistry,
coagulation & urine

*presets based on CLSI guidelines and tube manufacturer IFU

Rotor and Tube Holders Included

Features

- **Fixed Angle Rotor** - separate specimens at a fixed angle of 45 degrees
- **Easy To Use** - simple two-button operation to select a preset cycle or lock to a single cycle
- **6 Flex FA is Fully Programmable** - save & recall 10 easily customisable cycle settings
- **Lid Lights** - flash at the end of each cycle & continue flashing until the lid is opened
- **Automatic Lid Lock** - will not run with lid open & lid cannot be opened while samples are spinning
- **Cool Operation** - unique design of the convection based cooling system protects the samples even under the heaviest workloads
- **Clear Shatterproof Lid** - for easy monitoring
- **Brushless DC Motor** - no routine maintenance
- **Reinforced Guard Bowl** - maximum safety & durability

Specifications

	HORIZON 6 Flex FA	HORIZON 6 FA
Tube Capacity	6 tubes – 3-15 mL	6 tubes – 3-15 mL
Dimensions (H x W x D)	36 x 30 x 23 cm (14" x 12" x 9")	36 x 30 x 23 cm (14" x 12" x 9")
Weight	5.4kg (12 lbs)	5.4kg (12 lbs)
Sound Level	64 dB A	64 dB A
Environmental Range	16-32°C	16-32°C
Voltage	95-253 VAC	95-253 VAC
Frequency	50/60 Hz	50/60 Hz
Power Requirement	220 Watts	220 Watts
Centrifuge Motor	½ H.P. Brushless	½ H.P. Brushless
Max G-force	1,600 xg / 1,850 xg	1,600 xg / 1,850 xg
Max Speed	3900 RPM	3900 RPM
Cycle Time	1-30 minutes (+/- 2%)	1-30 minutes (+/- 2%)

Visit www.woodleylabdiagnostics.com for full product range